## **BENSAN PREMIUM SPRAY**



Safecare's Bensan Premium Spray mid-level alcohol solution is a ready-to-use spray for use on hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin. The formulation can be applied on polyethylene, propylene, stainless steel, plastic, Formica® and glass among others.

## **Efficacy Claims**

<ul> <li>Fungi/ Yeast incl. C. Auris</li> </ul>	1 Min.
Bacteria	1 Min.
Enveloped Viruses	1 Min.
<ul> <li>Non-Enveloped Viruses (incl. Noro)</li> </ul>	1 Min.
Mycobacterium Tuberculosis	1 Min.



- Compatible with a wide range of surface materials
- Free of phenols, bleach or aldehydes
- Available in 1-Liter bottles and 4-Liter canisters
- Effective cleaner with residual efficacy
- Non-toxic formulation
- Suitable for use on uneven hard-to-reach, porous and uneven surface areas
- Highly efficacious formulation based on Ethanol, Isopropyl Alcohol, 5th generation Quaternary Ammonium Compounds, ionic and non-ionic cleaning agents,
- Developed and manufactured in the United Arab Emirates with highest quality ingredient

## **List of Uses**

- Ambulance equipment surfaces
- Animal care facilities
- Correctional facilities
- Dental offices
- Emergency medical settings
- Emergency vehicles
- Exterior surfaces of anaesthesia machines and respiratory therapy equipment
- Hospitals
- Infant/child care equipment surfaces
- Isolation areas
- Laboratories
- Nursing homes
- Operating rooms
- Ophthalmic and optometric facilities
- Outpatient surgical centers
- Surgical centers

Bensan Premium in the same formulation is also available as a ready-to-use wipe.

ITEMS	REORDER NUMBER	CONTENT
Bensan Premium Spray	ICP2001	1 Liter
Bensan Premium Spray	ICP2002	4 Liters



Scan the QR Code for More Information

