



# Bensen Omni Concentrate

Safecare`s Bensen Omni alcohol-free concentrate is for the cleaning and disinfection of medical and dental instruments including rotating equipment prior to high-level disinfection or sterilization as well as the disinfection of environmental surfaces and surfaces of non-invasive medical devices.

Bensen Omni is highly suitable for instrument and equipment that are sensitive to acids, alkalis or alcohols. The concentrate consists of a combination latest generation quaternary ammonium compounds and surfactants and offers excellent cleaning and disinfection efficacy against a range of pathogenic microorganisms such as gram-positive and gram-negative bacteria, fungi and enveloped viruses within a 5-minute exposure time.

## Efficacy

- Fungi/ Yeast.....5 Min
- Bacteria:.....5 Min
- Enveloped Viruses.....3 Min
- Non-Enveloped Viruses.....5 Min

## Features

- Broad Spectrum
- Excellent material compatibility with sensitive medical instruments
- Excellent cleaner with bacteriostatic residual efficacy
- Strong concentrate, available in economic 4-Liter canisters – highly cost-effective

- Suitable for use on uneven hard-to-reach, porous and uneven surface areas
- Free of alcohol, bleach, phenols or aldehydes
- Non-toxic formulation
- Fragrance-Free or Fragranced (Lemon/Lavender) product available
- Biodegradable
- Follows US CDC instrument reprocessing guidelines
- Developed and manufactured in the United Arab Emirates with highest quality ingredients

### Uses:

- Hospitals
- Ambulance equipment surfaces
- Animal care facilities
- Correctional facilities
- Dental offices
- Emergency medical settings
- Emergency vehicles
- Exterior surfaces of anaesthesia machines and respiratory therapy equipment
- Infant/child care equipment surfaces
- Interior and exterior surfaces of infant incubators, bassinets
- Isolation areas
- Laboratories
- Nursing homes
- Operating rooms
- Ophthalmic and optometric facilities
- Outpatient surgical centers
- Surgical centers Bensen Omni in the same formulation is also available as a ready-to-use wipe.



Safecare Medical Industries  
 KIZAD KHIA 8-18  
 PO. Box 133685, Abu Dhabi  
 United Arab Emirates

For product or technical information, contact Safecare Medical at +971 2 506 7333 or visit our website at [www.bensano.com](http://www.bensano.com)