

Technical
Datasheet
10.08.2023

Bensan Omni Concentrate

Non-Alcoholic Disinfection Concentrate

*Instrument
Disinfection
Concentrate*



Application Areas

Safecare's Bensan Omni alcohol-free concentrate is for the cleaning and disinfection of medical and dental instruments including rotating equipment prior to high-level disinfection or sterilization.

Bensan Omni is highly suitable for instruments and equipment that are sensitive to acids, alkalis or alcohols.

Bensan Omni consists of a combination latest generation quaternary ammonium compounds and surfactants and offers excellent cleaning and disinfection efficacy against a range of pathogenic microorganisms such as gram-positive and gram-negative bacteria, and viruses.

Features:

- Strong concentrate – highly cost-effective (4 Liters yield up to 400 litres of Solution)
- Quick antimicrobial action – Virucidal (enveloped viruses), bactericidal, and fungicidal in 5 minutes
- Excellent cleaning function
- Excellent material compatibility with sensitive medical instruments
- Follows US CDC instrument reprocessing guidelines
- Fragrance free solution
- Biodegradable
- 3-in-1 cleaner, decontaminant and disinfectant

Composition

- 9.1% Various 5th Generation Quaternary Ammonium Compounds
- Surfactants
- Excipients

Shelf life

3 Years from date of manufacture

Efficacy Claims

Efficacy	Concentration	Contact Time	EN Norm
Virucidal	1% Dilution	5 Minutes	EN 14476
Bactericidal	1% Dilution	5 Minutes	EN 13727
Fungicidal	1% Dilution	5 Minutes	EN 13624, EN 13624

Chemical-physical data

Color	Colorless
Odor	Solvent odor
Flammability	Non-Flammable
Form	Liquid
pH	8-9

Toxicological Information

Carcinogenicity	None of the components of this product are listed as a carcinogen or suspected carcinogen by OSHA, IARC, and NTP
Chronic Effects	None known
Reproductive Effects	Not expected from this product
Mutagenic effects	Not expected to cause mutagenic activity
Eye	Can cause mild eye irritation
Skin	Non-irritating
Inhalation	High concentrations may cause upper respiratory tract irritation
Ingestion	No adverse effects are expected.

Best Practices

Use disinfectants safely. Always read the label and product information before use. Wearing protective personal equipment (e.g. made of nitrile rubber) is required. Do not mix with other cleaning agents. Further information is available on request.

Order Information

Name	Unit of Measurement (UOM)	Reorder Number
Bensan Omni (Concentrate)	10 X 1 L / case	ICO1006
Bensan Omni (Concentrate)	2 X 4 L / case	ICO1005

Related Products

- Bensan Omni Wipes
- Bensan Omni Spray

Safecare Medical Industries manufactures economically attractive high-quality disinfectants and instrument reprocessing products following environmentally friendly manufacturing processes in our state-of-the-art facility.

Safecare Medical Industries

KIZAD KHIA 8-18
PO. Box 133685, Abu Dhabi
United Arab Emirates
Telephone: +971 2 506 7333

