## Bensan Hypo Econo Wipes

Bensan Hypo Econo Wipes are presoaked ready-to-use intermediate level, sporicidal surface-disinfection wipes for use on hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin.

The Wipes can be applied on stainless steel, plastic, Formica® and glass among others.

The formulation consists of a combination of Hydrogen Peroxide, colloidal silver and surfactants and offers excellent residue-free

cleaning and disinfection efficacy against a wide range of pathogenic microorganisms within a 3 minute exposure time (5 minutes for spores). The non-toxic formulation which is compatible with a wide variety of medical devices can be applied broadly.

BENSAN

ensan Hypo Econo Wipe

mber: ICW1004

## **Examples include:**

- General patient areas
- Endoscopy
- Dialysis centers
- Operating rooms or
- CSSD

## Features:

- Sporicidal in 5 Minutes
- Completely residue-free
- Good material compatibility
- Completely non-toxic, biodegradable, environmentally friendly
- Effectively removes mold and bad odors
- Heavily soaked wipes

Efficacy	Contact Time	EN Norm
Virucidal (enveloped)	1 Min.	EN 14476
Virucidal (non-enveloped)	3 Min.	EN 14476
Bactericidal	3 Min.	EN 13727, EN16615
Fungicidal	2 Min.	EN 13624, EN16615
Tuberculocidal	3 Min.	EN 14348
Sporicidal	5 Min.	EN 17126

## **Order Information:**

Bensan Hypo Econo Wipes	10 canisters / case (80 wipes each)	ICW1004	

Contact info@bensano.com for a free trial or visit our website at www.bensano.com to request a free sample.

Safecare Medical Industries KIZAD KHIA 8-18 PO. Box 133685, Abu Dhabi United Arab Emirates Telephone: +971 2 506 7333 www.bensano.com