

Technical
Datasheet
15.02.2024
Rev. 1.0

Bensan Antifog

Antifogging Agent for Instrument Lenses, Goggles and Glasses

Bensan Antifog
Solution for lensed
instruments and
goggles



Application Areas

Bensan Antifog is intended for use during endoscopic, laparoscopic, gastroscopic, and arthroscopic procedures as well as any other procedures which require the use of an endoscope device, to prevent fogging of the endoscope lens. Packaged sterile. Single use.

Features:

- Designed to help prevent the scope or camera lens from fogging during laparoscopic, gastroscopic and arthroscopic surgeries
- Can also be used on protective eyewear, goggles and glasses
- Nontoxic and laboratory tested to be safe for use within a body cavity
- Includes 6 cc bottle with screw cap and adhesive-backed foam pad with x-ray detectable or radiopaque strip
- Not made with natural rubber latex
- Single use

Composition

- Water, Isopropyl Alcohol, Glycerin, Sodium Cumenesulfonate, Sodium lauryl Ether Sulfate, Citric Acid

Delivery Form

Form	Concentration
100 X 10 mL bottles	Ready-to-Use

Chemical-physical data

Color	Clear
Odour	Slight alcoholic odor
Storage	At room temperature in a dry, well ventilated area. Keep away from incompatible materials
Form	Liquid
pH Diluted Solution	7-8

Toxicological Information

Carcinogenicity	None of the components of this product are listed as a carcinogen or suspected carcinogen by OSHA and IARC.
Chronic Effects	None known

Reproductive Effects	None known
Mutagenic effects	Not expected to cause mutagenic activity
Eye	Rinse with plenty of water, get medical attention
Skin	None known. Flush with plenty of water during 15 minutes. Get medical attention.
Inhalation	None known (Remove to fresh air, get medical attention in case of symptoms)
Ingestion	Drink plenty of water. Get medical attention.

Material Compatibility

Stainless steel, anodized aluminum, chromed brass, nickel plated brass, ABS, polycarbonate, polypropylene, polystyrene and red natural rubber.

Shelf life

3 years from the date of manufacture

Instructions for Use

1. Dispense pouch contents onto sterile field.
2. Remove paper backing from foam pad and adhere foam pad onto sterile drape in a convenient location.
3. Remove cap and apply 5 or 6 drops of Anti-Fog solution to sponge.
4. Touch distal end of scope to sponge, blot dry (do not wipe dry) with sterile sponge.
5. Alternatively, apply solution directly to laparoscope lens.
6. Repeat as needed.

Order Information

Name	Unit of Measurement (UOM)	Reorder Number
Bensan Antifog	100 X 10mL Bottles	AVG1001

Safecare Medical Industries manufactures economically attractive high-quality disinfectants and instrument reprocessing products following environmentally friendly manufacturing processes in our state-of-the art facility.



Safecare Medical Industries

KIZAD KHIA 8-18
 PO. Box 133685, Abu Dhabi
 United Arab Emirates
 Telephone: +971 2 506 7333