

Value Spray by Safecare

Low-Alcohol Disinfection Spray

Surface Disinfectant Spray



Application Areas

Safecare's Bensen Value Spray low-alcohol solution is a ready-to-use spray for use on hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin. The formulation can be applied on stainless steel, plastic, Formica® and glass among others. The low-alcohol formulation combined with latest generation quaternary ammonium compounds offers excellent cleaning efficacy, allowing for a thorough disinfection in one pass within a 2 minute contact time. The non-toxic formulation which is compatible with a wide variety of medical devices and can be applied broadly. Examples include general patient areas, endoscopy, dialysis centers, operating rooms, CSSD or neonatal wards.

The spray can be used on any equipment that allows for a spray-and-wipe disinfection and is compatible with low levels of alcohol. Example include TEE probes, incubators, ultrasound probes, dialysis units, patient beds, mattresses, patient monitors among others.

The product is free of aldehydes, phenols or bleach.

Composition

- 17.6% Isopropanol
- 0.28% Quaternary Ammonium Compounds
- Surfactants

Packaging Sizes

Metric	1-Liter bottle
	4-Liter canister

Shelf Life

3 Years from the date of manufacture

Efficacy Claims

Efficacy	Concentration	Contact Time
Virucidal	Ready-to-use	2 Min. (limited virucidal)
Bactericidal	Ready-to-use	2 Min.
Fungicidal	Ready-to-use	2 Min.
Tuberculocidal	Ready-to-use	2 Min.

Chemical-physical data

Color	Colorless
Density	Approx. 0.97g / cm ³ at 20°C
Odor	Mildly alcoholic odor (isopropanol)
Flash Point	64°C
Form	Liquid
pH	10.5-12.5

Toxicological Information

Carcinogenicity	None of the components of this product are listed as a carcinogen
-----------------	---

	or suspected carcinogen by OSHA, IARC, and NTP
Chronic Effects	None known
Reproductive Effects	Not expected from this product
Mutagenic effects	Not expected to cause mutagenic activity
Eye	Can cause eye irritation
Skin	No adverse effects are known
Ingestion	No adverse effects are expected

Best Practices

Use disinfectants safely. Always read the label and product information before use. Wearing protective personal equipment (e.g. made of nitrile rubber) is required. Do not mix with other cleaning agents. Further information is available on request.

Order Information

Name	Unit of Measurement (UOM)	Reorder Number
Bensan Value Spray	10 X 1L / case	ICV1001
Bensan Value Spray	2 X 4L /case	ICV1002

Related Products

- Bensan Value Wipes
- Bensan Omni Spray and Wipes
- Bensan Premium Spray and Wipes
- Bensan Heavy Duty Spray and Wipes

Safecare Medical Industries manufactures economically attractive high-quality disinfectants and instrument reprocessing products following environmentally friendly manufacturing processes in our state-of-the art facility.

Safecare Medical Industries

KIZAD KHIA 8-18
PO. Box 133685, Abu Dhabi
United Arab Emirates
Telephone: +971 2 506 7333

