

Bensan IPA

by Safecare Medical Industries



Bensan IPA Surface is the latest addition to the extensive line of Bensan Surface Disinfectant Sprays, providing a 1-minute contact time from a 70% Isopropanol-based solution that fully complies with UAE and international standards.

Bensan IPA is free of quaternary ammonium compounds and provides a complete residue-free, non-allergenic solution to surface disinfection. The 70% alcohol are complemented by cleaning agents that provide a comprehensive 1-step – 1 minute reprocessing functionality for a broad spectrum of claims, effectively preventing protein fixation and enabling a deep and fast intermediate level disinfection.

USES

Bensan IPA should be used on all surfaces resistant to high-levels of alcohol, requiring a fast intermediate-level disinfection.

Bensan IPA can be used as a preliminary disinfection solution for surgical medical instruments and dental instruments prior to manual reprocessing to avoid cross-infection risks from needle-stick injuries and sharp cuts.

Bensan IPA can be used on areas hard to reach with a disinfectant wipe or on surfaces not suitable for use of a wipe. Make sure you soak surface areas sufficiently to allow for a 1-minute exposure time.

Bensan IPA has a 36-month shelf-life making it easy for bulk ordering and storage and comes in economical 4 liter-refill bottles.

EFFICACY CLAIMS

● The high-performance protocol of Bensan IPA offers a very broad highly efficacious solution against:

● **Enveloped viruses incl.**

- Herpes Simplex, Type 1 & Type 2
- Human Immunodeficiency Virus Type 1 (HIV-1)
- Hepatitis B virus (HBV)
- Hepatitis C Virus (HCV)

● Non-Enveloped viruses incl.

- Polio
- Human Norovirus
- Human Rotavirus

● Bacteria

- Pseudomonas Aeruginosa
- Staphylococcus Aurerus
- Salmonella Enterica

● Mycobacterium TB

- Mycobacterium Bovis (BCG)

● Multidrug-resistant Bacteria (MDRO)

- Penicillin-Resistant Streptococcus pneumoniae (PRSP)
- Vancomycin-Resistant Enterococcus faecalis (VRE)
- Methicillin-Resistant Staphylococcus Aureus (MRSA)

● Fungi

- Candida Auris
- Candida Albicans
- Trichophyton Mentagrophytes

● List of Uses

- Ambulance equipment surfaces
- Animal care facilities
- Bathrooms
- Correctional facilities
- Daycare centers
- Bordetella Pertussis
- Klebsiella Pneumoniae

Ordering Information

| Part No | Description | UOM |
|---------|---------------------------------------|---------|
| ICI1001 | Bensan IPA Surface Disinfectant Spray | 15 X 1L |
| ICI1004 | Bensan IPA Surface Disinfectant Spray | 4 X 4L |

Bensan IPA is entirely developed, produced and packaged in Safecare's factory in Abu Dhabi, UAE.

Safecare Medical Industries

KIZAD KHIA 8-18
PO. Box 133685, Abu Dhabi
United Arab Emirates
Telephone: +971 55 433 4560



To learn more about Bensan IPA, scan the QR code or visit us at www.bensano.com

