

Bensan Quinzyme

Quintuple-Enzymatic Detergent Concentrate for Manual and Automatic Use

Safecare's Bensan Quinzyme is an enzymatic detergent concentrate used for the manual and automatic cleaning of medical instruments. The 5 enzymes ensure a broad cleaning spectrum, leaving the instruments in crisp condition prior to sterilization or high-level disinfection. The formulation allows for short washing cycles, ensuring maximum throughput, faster turn-around time of critical instruments and reduced operational cost.

The formulation is based on protease, lipase, amylase, cellulase and mannanase as well as tensioactive agents. Thanks to the high enzymatic concentration and stable formulation, Bensan Quinyzme maintains a high level of enzymatic activity up to 24 hours after initial dilution providing further operational safety.

Features:

- Offering over 20% enzymatic activity, Bensan Quinzyme is one of the most highly concentrated enzymatic products available
- The dilution ranges from 0.05% to 0.5% providing up to 2000 liters of ready-to-use solution for each liter of Bensan Quinzyme
- Bensan Quinzyme offers excellent value and cost effectiveness
- Bensan Quinzyme consists of 5 enzymes, non-ionic surfactants, cleaning agents effectively removing bio-burden and biofilm build-up as well as calcium residue
- Bensan Quinzyme is pH neutral and includes corrosion inhibitors, providing excellent protection for the instruments
- The low-foaming formulation of Bensan Quinzyme offers ease of use during manual cleaning, using ultrasonic cleaners as well as in washer disinfectors and automatic endoscope reprocessors (AERs)
- The added sodium citrate and amine oxides offer strong anti-microbial activity maintaining bacteriostatic properties on gram-positive and gram-negative bacteria
- Bensan Quinzyme is biodegradable
- Bensan Quinzyme removes film build-up and residue left by high-level disinfectants
- On request and on meeting minimum purchase requirements, Safecare will provide free washer disinfector efficacy test strips

Description	Article Number	Content
Bensan Quinzyme 4L	IRQ1001	2 X 4 Liters/Case
Bensan Quinzyme 1L	IRQ1002	10 X 1 Liter/Case
Bensan Quinzyme 10L	IRQ1005	1 X 10 Liters
Bensan Quinzyme 20L	IRQ1020	1 X 20 Liters







For product or technical information, contact



