

# Bensan Bleach Econo Wipes

Bensan Bleach Econo Wipes are presoaked, ready-to-use wipes for use on hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin. The formulation can be applied on stainless steel, plastic, Formica® and glass among others.

The bleach wipes offer excellent cleaning efficacy, allowing for a thorough disinfection in one pass within a 3 minute contact time (1 minute for enveloped viruses and bacteria). The wipes which are compatible with a wide variety of medical devices and can be applied broadly. Examples include general patient areas, endoscopy, dialysis centers, operating rooms, CSSD or neonatal wards.

Bensan Bleach Econo Wipes are suitable for use on any surface material that is compatible with bleach.

## Example include:

- TEE probes
- Incubators
- Ultrasound probes
- Dialysis units
- Patient beds
- Mattresses
- Patient monitors among others
- Counter tops
- Keyboards

## Features:

Bactericidal, Virucidal, Fungicidal and Tuberculocidal
Phenol-free, aldehyde-free formulation
Low odor
Non-toxic
Highly saturated wipes
Large coverage area / wipe – highly cost-effective
80 Wipes of 30 gsm / canister
Developed and Manufactured 100% in the UAE

Efficacy	Concentration	ontact Time	EN Norm
Virucidal (enveloped)	Ready-to-use	3 Min.	EN 14476
Virucidal (non-enveloped)	Ready-to-use	3 Min.	EN 14476
Bactericidal	Ready-to-use	1 Min.	EN 13727, EN16615
Fungicidal	Ready-to-use	1 Min.	EN 13624, EN16615
Tuberculocidal	Ready-to-use	3 Min.	EN14348

## Order Information:

Bensan Bleach Econo Wipes	Bensan Bleach Econo Wipes	BBW1080
---------------------------	---------------------------	---------

Contact [info@bensano.com](mailto:info@bensano.com) for a free trial or visit our website at [www.bensano.com](http://www.bensano.com) to request a free sample.

Safecare Medical Industries  
 KIZAD KHIA 8-18  
 PO. Box 133685, Abu Dhabi  
 United Arab Emirates  
 Telephone: +971 2 506 7333  
[www.bensano.com](http://www.bensano.com)

