

Technical  
Datasheet  
03.11.2023  
Rev. 2.0

# Bensan Premium Econo Wipes Alcoholic Disinfection Wipes

*Alcoholic Surface  
Disinfectant  
Wipes*



## Application Areas

Bensan Premium Econo Wipes mid-level alcoholic presoaked wipes are for use on hard, non-porous environmental surfaces as well as surfaces of non-critical medical devices that get in touch with the intact skin. The wipes can be applied on stainless steel, plastic, Formica®, Polyethylene, Polypropylene, Polyvinyl Chloride and glass among others. The formulation consists of a combination of Isopropanol, Ethanol and a latest generation quaternary ammonium compounds offers excellent cleaning efficacy and disinfection efficacy against a wide range of pathogenic microorganisms within a 1 minute exposure time. The non-toxic formulation which is compatible with a wide variety of medical devices can be applied broadly. Examples include general patient areas, endoscopy, dialysis centers, operating rooms, CSSD or neonatal departments.

The solution can be used on any equipment that allows for a spray-and-wipe disinfection and is compatible with medium levels of alcohol. Example include TEE probes, incubators, ultrasound probes, dialysis units, patient beds, mattresses, patient monitors among others.

All products developed and manufactured in the United Arab Emirates.

## Composition

• 27.3% Isopropanol
• 28.7% Ethanol
• 0.61% Quaternary Ammonium Compounds
• Surfactants, Excipients, corrosion inhibitors, inert ingredients

## Packaging Sizes

• 10 canisters per case	80 Wipes/Canister
-------------------------	-------------------

## Shelf Life

3 Years from the date of manufacture

## Efficacy Claims

Efficacy	Concentration	Contact Time	EN Norm
Virucidal	Ready-to-use	1 Min.	EN 14476
Bactericidal	Ready-to-use	1 Min.	EN 13727, EN 16615
Fungicidal	Ready-to-use	1 Min.	EN 13624, EN 16615
Yeasticidal	Ready-to-use	1 Min.	EN 13624, EN 16615
Tuberculocidal	Ready-to-use	1 Min.	EN 14348

## Chemical-physical data

Color	Colorless
Density	Approx. 0.97g / cm <sup>3</sup> at 20°C
Odor	Alcoholic odor
Flammable	Keep away from heat or open flame
Form	Liquid

pH	8.5- 11.5
----	-----------

## Toxicological Information

Carcinogenicity	None of the components of this product are listed as a carcinogen or suspected carcinogen by OSHA, IARC, and NTP
Chronic Effects	None known
Reproductive Effects	Not expected from this product
Mutagenic effects	Not expected to cause mutagenic activity
Eye	Can cause eye irritation
Skin	Prolonged exposure can cause moderate irritation with drying and dermatitis
Inhalation	High concentrations may cause upper respiratory tract irritation
Ingestion	Ingestion is unlikely for solid products. No adverse effects are expected.

## Best Practices

Use disinfectants safely. Always read the label and product information before use. Wearing protective personal equipment (e.g. made of nitrile rubber) is required. Do not mix with other cleaning agents. Further information is available on request.

## Order Information

Name	Unit of Measurement (UOM)	Reorder Number
Bensan Premium Econo Wipes	20 canisters / case (80 wipes each)	ICP2003

## Related Products

- Bensan Value Econo Wipes
- Bensan Omni Spray and Omni Econo Wipes
- Bensan Premium Spray
- Bensan Bleach Econo Wipes
- Bensan IPA Econo Wipes

**Safecare Medical Industries** manufactures economically attractive high-quality disinfectants and instrument reprocessing products following environmentally friendly manufacturing processes in our state-of-the-art facility.

## Safecare Medical Industries

KIZAD KHIA 8-18  
PO. Box 133685, Abu Dhabi  
United Arab Emirates  
Telephone: +971 55 433 4560

