

# Bensan iTreat

## Povidone Iodine-based Antiseptic Tincture

*Povidone Iodine &  
Ethanol-based  
Antiseptic Tincture*



### Application Areas

Safecare's Bensan iWash antimicrobial wound treatment solution is based on a highly efficacious yet skin-friendly formulation of povidone iodine and ethanol. The solution is highly effective against viruses (enveloped and non-enveloped), bacteria, fungi and TB. The residual anti-microbial properties hours make the solution ideal for the strong demands of a wound treatment tincture as well as skin prep for surgical use. The solution does not interact with Latex, Vinyl or Nitrile Gloves and allows for quick donning after use. Emollients maintain skin in healthy condition.

The solution complies with the EU Biocides Regulation 528/2012.

### Composition

- 50% Ethanol
- 10% Povidone Iodine (1% active Povidone Iodine)
- Sodium Iodide
- Water

### Shelf life

3 Years from the manufacturing date

### Efficacy Claims

Efficacy	Concentration	EN Norm
Virucidal – Enveloped Viruses	Ready-to-use	prEN 1476, EN1500
Topical Bactericide	Ready-to-use	prEN 12054
Fungicidal	Ready-to-use	EN 1500

## Chemical-physical data

Color	Dark brown
Odor	Iodine odor
Flammable	No
Form	Liquid
pH	5-7

## Toxicological Information

Carcinogenicity	None of the components of this product are listed as a carcinogen or suspected carcinogen by OSHA, IARC, and NTP
Chronic Effects	None known
Reproductive Effects	Not expected from this product
Eye	Can cause eye irritation
Skin	No adverse effects are known
Inhalation	No adverse effects are known
Ingestion	No adverse effects are expected



## Instructions for Use

- Gently clean the affected area with cotton soaked in the solution and 2.5 cm of the peri-wound skin
- Allow the solution to dry
- Cover the wound with a sterile bandage
- Bottle must be client-single use, labeled with client's name and date of bottle

## Hazard Statements



Causes Serious Eye Irritation



Highly Flammable liquor or vapour

## Order Information

Name	Unit of Measurement (UOM)	Reorder Number
Bensan iTreat	20 X 500 mL / case	PVP2001

**Safecare Medical Industries** manufactures economically attractive high-quality disinfectants and instrument reprocessing products following environmentally friendly manufacturing processes in our state-of-the-art facility.

## Safecare Medical Industries

KIZAD KHIA 8-18  
 PO. Box 133685, Abu Dhabi  
 United Arab Emirates  
 Telephone: +971 2 506 7333  
[www.bensano.com](http://www.bensano.com)  
[info@bensano.com](mailto:info@bensano.com)

