



# SafeWipe Dry



## Quintuple-Action

- Precleaning using a dry towel
- Cleaning with a medical grade detergent (spray-and-wipe)
- Decontamination with a low- to intermediate level disinfection (spray-and-wipe)
- Disinfection with an intermediate level disinfectant (spray-and-wipe)
- Lint-free and aseptic drying pass all with the same product

## FEATURES

Allow for a complete reprocessing of a wide range of surfaces

Effectively remove soil and bacteria

Large, dense and abundant wipes 32X30 cm

Lint-free, residue-free

Ubiquitous application

Highly absorbent

Ideal for use with all surface reprocessing chemicals

Ideal for professional and personal use

Economical packs of 50 towels/pack

SafeWipe Dry quintuple-action towels are made of highly absorbent lint-free microfibers specifically configured and manufactured for the optimized reprocessing of environmental surfaces and surfaces of medical equipment.

Safecare's maximum-efficacy SafeWipe Dry towels eliminate mechanically 99.998% of bacterial contamination preparing surfaces for a comprehensive, thorough reprocessing cycle.

When combined with suitable detergent and disinfectant, SafeWipe Dry towels offer a comprehensive surface reprocessing solution:

SafeWipe Dry towels' abundant, large size of 32cm by 30cm provides ample surface coverage of up to 2 square meters making them highly cost effective.

Protect yourself, your family, your friends, your staff and your patients with Safecare Medical's SafeWipe Dry towels.



**Safecare Medical Industries**  
KIZAD KHIA 8-18, PO. Box 133685,  
Abu Dhabi United Arab Emirates

**Item:** SafeWipe Dry

**Reorder Number:** CSAS1005

**UOM:** 10 X 50 Towels/Case